



Products / Traditional Wound Care / Non-Woven Dressings / Non-Woven Sponges /  
Non-Sterile Non-Woven Sponges, 4-Ply, 2" x 2"



## Non-Sterile Non-Woven Sponges, 4-Ply, 2" x 2"

2102

200/BG, 20 BG/CS

### Product Features

- Rayon/poly blend construction
- Highly absorbent
- Multi-purpose sponges ideal for general wound care and applying ointments
- Virtually non-linting

### Attributes

Brand Name	Dukal
Size	2" x 2"
Ply	4-Ply
HCPC Code	A6216
Sterile	No
Weight	2.95 lbs
Item Status	STOCK

### Packaging

Packaging	200/BG, 20 BG/CS
Inner Packaging	200/BG
Case Packaging	20 BG/CS
Total Case Quantity	4000
Inner Packaging Barcode	(01)10665793013760
Case Packaging Barcode	(01)20665973013767

2 Fleetwood Ct.  
Ronkonkoma, NY 11779

Phone: 1 (631) 656-3800  
Toll Free: 1 (800) 243-0741  
Fax: 1 (631) 656-3810

