

Modular Special Procedures Cart with Two Drawers

Model 6620

The UMF Medical 6620 special procedures cart has two drawers, a stainless steel push handle, and 5-leg base with twin wheel carpet casters for maximum stability. It comes equipped with a Kydex work surface and is the perfect complement to any exam room or treatment area.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel construction
- › Stainless steel push handle
- › Kydex top with rail
- › 2" (5.1cm) twin wheel carpet casters
- › 3-year limited warranty

Specifications

Width: 20" (50.8cm)
Height: 33.5" (85.1cm)
Depth: 16.25" (41.3cm)
Net Weight: 49lb (22.2kg)

Shipping Weight: 62lb (28.1kg)
Shipping Dimensions:
L21xW28xH22 (53.3cm x
71.1cm x 55.9cm)

Accessories

116 – Swivel casters locking 60mm
(factory installed) (set of 5)
505 – Single drawer lock (factory installed)
507 – Additional keys (x2)