



Products / Surgical Products / Non-Woven Dressings / T-Drain Sponge / Sterile Non-Woven T-Drain Sponges 4" x 4" 6-Ply



Sterile Non-Woven T-Drain Sponges 4" x 4" 6-Ply

7470

2/PK, 25 PK/BX, 12 BX/CS

Product Features

Pre-cut rayon/poly blend sponges

Designed to wick moisture away from the site

Unique pre-cut T-slit provides a snug fit around catheters, tubes, I.V.'s and tracheotomies

T-Drain sponge also available in 100% woven cotton

Attributes

Brand Name	Dukal
Size	4" x 4"
Ply	6-Ply
HCPC Code	A6402
Sterile	Yes
Weight	5.75 lbs
Item Status	STOCK

Packaging

Packaging	2/PK, 25 PK/BX, 12 BX/CS
Item Packaging	2/PK
Inner Packaging	25 PK/BX
Case Packaging	12 BX/CS
Total Case Quantity	600
Item Packaging Barcode	(01)10665973016389
Inner Packaging Barcode	(01)30665973016383
Case Packaging Barcode	(01)20665973016386

2 Fleetwood Ct.
Ronkonkoma, NY 11779

Phone: 1 (631) 656-3800
Toll Free: 1 (800) 243-0741
Fax: 1 (631) 656-3810

