

# Bariatric & Parent with Child Phlebotomy Chair

## Model 8690

UMF Medical's 8690 Bariatric Blood Draw, "Mommy-and-Me" Chair offers a flip-up, height-adjustable arm board & front board. Heavy-duty steel frame, with a 750lb weight capacity. Available in all 16 upholstery colors with PreFixx® Protective Finish.



### Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



### Key features at a glance

- › Ultra comfort seamless upholstery
- › Accommodates bariatric patients and parent with child
- › Padded, flip-up, height adjustable arm board
- › Heavy duty, all-welded steel frame
- › 750lb (340.2kg) patient weight capacity
- › Adjustable leg levelers
- › 3-year limited warranty

### Specifications

Width: 41" (104cm)  
 Depth: 35" (88.9cm)  
 Height: 41.75" (106cm)  
 Net Weight: 78lb (35.4kg)  
 Shipping Weight: 138lb (62.6kg)  
 Shipping Dimensions: L44xW34xH50  
 (112cm x 86.4cm x 127cm)

Seat dimension:  
 Height: 22" (55.9cm)  
 Width: 26.625" (68cm)  
 Between arm boards:  
 30" (76cm)

### Accessories

120 – Side tray  
 121 – Storage Box