

THE PIONEER OF MEDICAL
EDUCATIONAL SIMULATOR +

Blood Pressure Assessment Simulator

BT-CEAB2

The various case of disease and exact Real-Time Feedback about the training

NURSING MEDICINE MILITARY EMS HOSPITAL



BP/ Pulse Assessment Simulator

MANIKIN

- Palpate pulses at Brachial/Radial Artery
- Group Training for Korotkoff Sound with Auxiliary Speakers
- Tablet PC ▶ Evaluate practical training



Pulse Sites



Radial



Brachial



Korotkoff Sounds



Brachial




Speaker

SOFTWARE

- Mode : Auscultation / Palpation / Auscultation-Palpation
- Auscultatory gap training
- Detailed feedback in real-time
- Training Mode / Evaluation Mode
- Objective evaluation with a manual checklist: Scoring

- Systolic BP
- Diastolic BP
- Pulse Rate
- Cuff Installation
- Palpation of Brachial Artery
- Auscultatory Gap
- Keeping time limit
- Proper compression
- Normal Decompression Rate

- Saving and Printing out results



* Display



* Setting



* Evaluation Mode Report

COMPONENTS

- Adult Arm Manikin : 1EA
- Tablet : 1EA
- Arm model mount : 1EA
- Cuff : 1EA
- Storage case : 1EA
- Power Adaptor : 1EA / Battery holder : 1EA



SPECIFICATION

- Adult arm model : 72 x 11 x 11cm (1.2kg)
- Display device : 26 x 18 x 1cm (0.6kg)
- Storage case : 73 x 28 x 14cm (2.2kg)