

THE PIONEER OF MEDICAL
EDUCATIONAL SIMULATOR +

Male/Female Catheterization Training Model

BT-CSCM / BT-CSCF

NURSING MEDICINE MILITARY EMS HOSPITAL



Male/Female Catheterization Training Model

FEATURES

- Life-like practical training
- Make position similar to the real patient
- Any water-soluble Jelly being used at clinical situation can be applied
- Training with various size of catheter is available
- Non-return valve : Prevent leaking after taking catheter out.
- Genital organs and buttocks are attachable : easy to maintain

Male - CSCM

- Foreskin (attachable / detachable)
- Checking fluid after correct insertion of catheter at the right position



With Foreskin

Without Foreskin



Female - CSCF

- Realistic practice with attachable extended legs
- Alarm when a catheter inserted into a vagina ; to operate, a battery is needed
- Checking fluid after inserting a catheter into bladder correctly



Warning Alarm



Extended legs

COMPONENTS

Male - CSCM

- Adult Buttocks model : 1 SET
- Simulated urine bag : 1EA
- Storage case : 1EA

Female - CSCF

- Adult Buttocks model : 1 SET
- Extension legs : 1 SET
- Simulated urine bag : 1EA
- Storage case : 1EA

SPECIFICATION

Male - CSCM

- Adult buttocks model : 3.5kg (430×280×200mm)
- Storage case : 2.4kg (440×330×390mm)
- Total : Approximately 6kg

Female - CSCF

- Adult buttocks model : 3.5kg (430×280×200mm)
- Extended legs (Left/Right) : Each 1kg (320×150×150mm)
- Storage case : 2.4kg (440×330×390mm)
- Total : Approximately 8kg