

THE PIONEER OF MEDICAL
EDUCATIONAL SIMULATOR +

Enema Administration Simulator

BT-CSEN

NURSING MEDICINE MILITARY EMS HOSPITAL



FEATURES

- Lifelike practical training
- Position similar to the real patient
 - Sim position / knee-chest position
- Real Glycerin can be used
- Real water-soluble jelly can be used
- Non-return valve : Prevent leaking after glycerin enema training
- Artificial feces can be inserted / removed



Using real Glycerin & Jelly



Enema Training



Artificial Feces

COMPONENTS

- Adult buttock model : 1EA
- Stand : 1EA
- Drain tube : 1EA
- Cleaning Syringe : 1EA
- Storage case : 1EA

SPECIFICATION

- Adult Buttocks model : 3.5kg (330×430×340mm)
- Storage case : 2.4kg (440×330×390mm)
- Total : Approximately 6kg