

THE PIONEER OF MEDICAL
EDUCATIONAL SIMULATOR +

Sherpa, Self-Training CPR Model

BT-SEEM-P

CPR Training with real-time check via abdomen LED

NURSING MEDICINE MILITARY EMS HOSPITAL



Sherpa - Self-training CPR model

FEATURES

- Chest compression depth and rate can be checked with LED lamp
- Training for Chest Compression
- Count total and correct number of compressions
 - Compression Depth and Speed
 - Instant compression depth
 - Incomplete recoil
 - Over-compression



*Display of Sherpa



*Compression Depth



*Compression Speed

- Sound and LED Lamp can be turned on/off for an advanced training
- Visual chest rising during ventilation
- Self replacement for face skin and lung
- Semi-permanent use of AED training pads with embedded magnets in manikin
(Offering magnetic stickers)



COMPONENTS

- Adult Torso : 1EA
- Storage case (use for mattress) : 1EA
- Power Adaptor : 1EA
- Magnetic stickers : 1SET (2EA)

SPECIFICATION

- Dimensions : 65 x 33 x 20cm
- Weight : 4.2kg

