



Shown with optional stainless steel wire shelves

CDS-150

Sterilizer Cart

Overall size: 37³/₄"W x 24¹/₂"D x 45³/₈"H
(includes bumpers)



Shown with optional stainless steel wire shelves

CDS-151

Sterilizer Cart

Overall size: 37³/₄"W x 24¹/₂"D x 55³/₈"H
(includes bumpers)

sterilizer carts

An excellent cart for autoclaves due to single piece, all welded stainless steel construction

features

- Fully-welded seams
- Heavy-duty 1¹/₄" diameter tubular frame
- 5" high temperature stainless steel swivel plate casters for smooth mobility
- High temperature corner bumpers
- Open wire bottom for improved steam circulation
- 17¹/₂" center-to-center caster spacing