



**CDS-256-LD**

Overall outside dimensions:  
28<sup>5</sup>/<sub>8</sub>"W x 28<sup>5</sup>/<sub>8</sub>"D x 29<sup>7</sup>/<sub>8</sub>"H

Usable Interior Space  
20<sup>1</sup>/<sub>4</sub>"W x 23<sup>1</sup>/<sub>2</sub>"D x 18<sup>1</sup>/<sub>2</sub>"H

*One shelf configuration recommended*



CDS-256-MS shown with optional drop leaf

**CDS-256-MS**

Overall outside dimensions:  
28<sup>5</sup>/<sub>8</sub>"W x 28<sup>5</sup>/<sub>8</sub>"D x 39"H

Usable Interior Space  
20<sup>1</sup>/<sub>4</sub>"W x 23<sup>1</sup>/<sub>2</sub>"D x 27<sup>1</sup>/<sub>2</sub>"H

*One shelf configuration recommended*

# case carts

surgical

**features**

- Fully-welded all stainless steel construction
- Exterior push handles for better control and more interior space
- One-piece continuous wrap-around bumper, non-marring
- Sound-deadened tops and doors for quiet operation
- Positive lock door latches
- All carts can be manufactured to individual specifications



Biohazard Indicator Wheel  
**CartStat**



Wire Shelf



6" Casters



**re•volve™**  
caster upgrade



CDS-235 shown with optional shelves & CartStat



CDS-242 shown with optional shelves



CDS-245 shown with optional shelves

**CDS-233**

Overall outside dimensions:  
28<sup>5</sup>/<sub>8</sub>"W x 28<sup>5</sup>/<sub>8</sub>"D x 56"H

Usable Interior Space  
20<sup>1</sup>/<sub>2</sub>"W x 23<sup>1</sup>/<sub>2</sub>"D x 44<sup>3</sup>/<sub>4</sub>"H

*Two shelf configuration recommended*

**CDS-235**

Overall outside dimensions:  
34<sup>1</sup>/<sub>8</sub>"W x 32<sup>3</sup>/<sub>8</sub>"D x 56"H

Usable Interior Space  
25<sup>3</sup>/<sub>4</sub>"W x 27<sup>3</sup>/<sub>8</sub>" x 44<sup>3</sup>/<sub>4</sub>"H

*Two shelf configuration recommended*

**CDS-242**

Overall outside dimensions:  
45<sup>7</sup>/<sub>8</sub>"W x 28<sup>5</sup>/<sub>8</sub>"D x 39"H

Usable Interior Space  
37<sup>1</sup>/<sub>2</sub>"W x 23<sup>1</sup>/<sub>2</sub>"D x 27<sup>3</sup>/<sub>4</sub>"H

*One shelf configuration recommended*

**CDS-245**

Overall outside dimensions:  
45<sup>7</sup>/<sub>8</sub>"W x 28<sup>5</sup>/<sub>8</sub>"D x 56"H

Usable Interior Space  
37<sup>1</sup>/<sub>2</sub>"W x 23<sup>1</sup>/<sub>2</sub>"D x 44<sup>3</sup>/<sub>4</sub>"H

*Two shelf configuration recommended*

**Shelving Configurations**

**CDS-233-SS** ..... Stainless steel solid shelf for CDS-233, CDS-256-LD and CDS-256-MS.

**CDS-233-SSRO** ... Stainless steel roll-out solid shelf for CDS-233, CDS-256-LD and CDS-256-MS.

**CDS-233-WS** ..... Stainless steel wire shelf for CDS-233, CDS-256-LD and CDS-256-MS.

**CDS-233-WSRO** .. Stainless steel roll-out wire shelf for CDS-233, CDS-256-LD and CDS-256-MS.

**CDS-235-SS** ..... Stainless steel solid shelf for CDS-235.

**CDS-235-SSRO** ... Stainless steel roll-out solid shelf for CDS-235.

**CDS-235-WS** ..... Stainless steel wire shelf for CDS-235.

**CDS-235-WSRO** .. Stainless steel roll-out wire shelf for CDS-235.

**CDS-242-SS** ..... Stainless steel solid shelf for CDS-242 and CDS-245.

**CDS-242-SSRO** ... Stainless steel roll-out solid shelf for CDS-242 and CDS-245.

**CDS-242-WS** ..... Stainless steel wire shelf for CDS-242 and CDS-245.

**CDS-242-WSRO** .. Stainless steel roll-out wire shelf for CDS-242 and CDS-245.

Custom CDS carts engineered to customer specifications available upon request.