

distribution carts



Cart shown with optional accessories

features

- All welded stainless steel construction
- All seams fully welded
- Frame tubing is 1 1/4" - 18 gauge stainless steel
- Wrap around bumper mounted away from the frame for easy cleaning
- Shelves adjust every four inches
- Low maintenance
- Solid stainless steel pan bottom with reinforcing channel and sound control strip for added strength and quiet operation

CDS-147 distribution Cart

OVERALL SIZE

39 3/8"W x 26 1/8"D x 42 1/8"H

- Ideal for case cart systems or general floor use
- Shelf size 35 1/4"W x 20 3/8"D
- 5" chrome plated swivel casters



CDS-147-A Swivel Caster



Cart shown with optional accessories

CDS-148 distribution Cart

OVERALL SIZE

50 5/8"W x 26 1/8"D x 58"H

- Suitable for general floor service, pharmacy or linen
- Shelf size 46 3/8"W x 20 3/8"D
- 8" chrome plated casters with two diagonally mounted swivel locks



Cart shown with optional accessories


CDS-149 distribution Cart

OVERALL SIZE

60 3/8"W x 26 1/8"D x 67"H

- Ideal cart for linen distribution, storage and floor nursing stations
- Shelf size 55 3/4"W x 20 3/8"D
- 8" chrome plated casters with two diagonally mounted swivel locks

Accessories



Retaining Rods (adjustable)

- Single style
- Double style



Tote Boxes (two way)

- 5" deep
- 10" deep




Shelf Dividers (adjustable)

- 6" tall
- 10" tall



Wire Backs

- 2" x 3" openings
- 4" x 6" openings



Stainless Steel Shelves

- Wire shelf
- Solid shelf