



### CDS-1100

Main Table Top Surface  
72"W x 30"D

Height Adjustable  
32" to 43¼"H

## instrument assembly station table

### features

- Electrically operated hydraulic height adjustment provides comfortable, ergonomic height settings for all staff members.
- Shelves are positioned 18" and 30" above main work surface, and measure 72"W x 12"D with 5/8" front lip.
- Shelves also feature a 10° slope.
- Load capacity of 225 pounds on main work surface and 50 pounds on each upper shelf.
- 4" zinc plated steel brake/swivel lock casters with precision ball bearings for easy mobility.

### accessories

#### 1 Work Lamp **CDS-1020**

- Work Lamp arm has a 42" reach
- Features a 5" viewing lens
- 75% magnification

#### 2 Outlet Strip **CDS-1015**

- Hospital grade
- Six outlets
- 15' foot cord

#### 3 Keyboard Tray **CDS-1032**

- Measures 22¼"W x 10"D x 2¾"H
- Tray extends 10½" (retractable)
- Built-in mouse tray comes standard

#### 4 Computer Enclosure **CDS-1031**

- Open back for air flow and easy access to computer cords.
- Enclosure can be customized to accommodate larger computer towers (max. load 30 lbs/13.60 kg)

#### 5 Extension Arm **CDS-1030**

- Supports LCD monitors (max. load 25 lbs / 11.3 kg)
- For optimal viewing, extension arm features multi-directional adjustment
- Gas-strut and ball-joint system for smooth, stable operation

#### 6 Storage Cart **CDS-1010**

- 4 drawer tote box cart for added storage fits under the workstation
- Totes are constructed of impact-resistant, high-density polyethylene and measure 17¾"L x 14½"W x 5"D