

DANi Teen & Young Adult Advanced Training Solution

Product Name & Overview	DANi Teen & Young Adult Advanced Training Solution
Addresses Training Needs/Curricula	<ul style="list-style-type: none">• Basic Life Support (BLS)• Advanced Cardiac Life Support (ACLS)• Certification in Sexual Abuse Trauma (CSAT)• Sexual Assault Nurse Examiner (SANE)• Sexual Assault Nurse Examiner-Pediatric/Adolescent (SANE-P)
Configuration	<ul style="list-style-type: none">• Configurable to fit unique training needs and budgets• A range of accessories available to maximize versatility
Consumables	<p>IV arm:</p> <ul style="list-style-type: none">• Skins• Veins kit• IM pads• Fluid supply bag• Simulated blood <p>Urinary catheterization:</p> <ul style="list-style-type: none">• Male genitalia foreskins• Foley catheter• Urinary catheterization bladders• Simulated urine <p>Additional consumables:</p> <ul style="list-style-type: none">• Lung bags• IM pads for shoulder, thighs and buttocks• Stomach/stoma/gastronomy bags• Simulated blood• Lubricant• Stomach bags• Enema reservoir
Product Code/SKUs	<p>Patient Care Trainer LF04500 (Light skin tone) LF04500M (Medium skin tone)</p> <p>Patient Skills Trainer LF04502 (Light skin tone) LF04502M (Medium skin tone)</p> <p>Auscultation Patient Skills Trainer LF04503 (Light skin tone) LF04503M (Medium skin tone)</p> <p>Patient Simulator LF04504 (Light skin tone) LF04504M (Medium skin tone)</p>
Dimensions/Weight	Dimensions: 22.25” H x 44.5” L x 25.5” W Ship Weight: 80 lbs.



Nasco Healthcare Inc.
16 Simulaids Drive, Saugerties, NY, 12477 USA
Phone: 1-833-NASCOHC (627-2642) Email: info@nascohealthcare.com
www.nascohealthcare.com



Be READY for the next generation of youth medical simulation.

Evolving from our popular KERi™ manikin, DANi is an advanced training solution simulating a teen/young adult patient with enhanced aesthetics for added realism.

DANi features superior range of motion, with weight distributed to accurately represent a real patient for lifting and carrying techniques.

A realistic, hands-on learning experience that is ideally suited for classroom, in-situ, or hospital instruction.

Each version is affordable and durable with upgradable items.

Patient Care Trainer

General patient care, patient positioning, bandaging, wound dressing, oxygen administration.

Patient Skills Trainer

Respiratory support, injection/IV with IM pads, intubation, CPR with quality performance metrics, non-invasive blood pressure, enema administration (female), catheterization, vaginal, pelvic and prostate examinations.

Auscultation Patient Skills Trainer

All the features of Patient Skills Trainer with added auscultation. Lung sounds at five anterior, ten posterior, and two midaxillary locations. Six anterior heart sites.

Patient Simulator

Comprehensive and complete simulation platform including patient monitoring capability and auscultation (heart and lung sounds).



DANi Teen & Young Adult Advanced Training Solution

Configurable to fit unique training needs and budgets with a range of accessories available to maximize versatility.

Trauma Simulation¹⁻⁴

- Physical assessment
- Airway, breathing, circulation
- Body trauma scan
- Wound care
- Moulage
- Asymmetric pupils




Realistic Features¹⁻⁴

- Pupils
- Airway
- Acne
- Smooth skin
- Tattoo sleeve
- Piercings & earrings
- Wigs
- Foot wounds

Patient Positioning¹⁻⁴

Superior range of motion, seating, laying, finger and toe manipulation



Nursing Care¹⁻⁴

Grooming, daily living assistance, eye/ear/hair care, catheter, injection, oxygen, blood pressure, pericare, ostomy care, urinary catheter care, oxygen administration

Respiratory Support²⁻⁴

Manual ventilation with chest rise, respiratory support, intubation, ventilation, tracheostomy care & suctioning

Auscultation^{3,4}

- 6 anterior heart sites with 12 heart conditions
- 5 anterior, 6 upper posterior, 4 lower posterior, 2 midaxillary lung sites with 12 lung conditions

Basic Life Support (BLS) & Advanced Cardiac Life Support (ACLS)²⁻⁴

CPR with quality performance metrics, carotid pulse, intubation

Advanced Skills²⁻⁴

Gastrostomy, nasogastric procedures, airway/trach suctioning, enema, pap smear, pelvic exam, intravenous injections

Female & Male Genitalia Inserts¹⁻⁴

REFERENCE GUIDE:

1. Available on Patient Care Trainer
2. Available on Patient Skills Trainer
3. Available on Auscultation Patient Skills Trainer
4. Available on Patient Simulator

This manikin is intended to be used as an instructional aide for teaching purposes. Actual product may vary slightly from photo. Nasco Healthcare reserves the right to change product color, materials, supplies, or function as needed.