



## Covered Square Hamper

- Square rolling hamper with easy-open lid
- 3" Black swivel casters: two locking, two non-locking
- Black metal footpad with rubber pad.
- Chrome-plated frame
- 16" x 12" hamper opening
- Optional Hamper Bag (not included with hamper)

<b>5532B</b>	Hamper	1 ea
<b>5532B-BAG</b>	Optional Hamper Bag	1 ea

Item	Height	Width	Depth
<b>5532B</b>	38.5"	19"	21"



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**