

The Advantage Wheelchair provides comfort and durability at an affordable price. Designed as a low maintenance wheelchair, it offers an attractive silver vein durable hammertone steel finish.



Standard Chart Pocket

Standard Features

- Available in 16", 18" and 20" seat widths
- Black padded nylon and vinyl upholstery available with standard chart pocket
- Carbon steel frame provides strength and durability
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Hammertone paint provides a long lasting and durable finish
- Comfortable padded armrests
- Maximum Weight Capacity 300 lb – EVENLY DISTRIBUTED
- Meets Medicare code K0001

Warranty

One (1) year limited warranty



**This product is not made with natural rubber latex*












Fixed Full Arm,
Swingaway Footrest



Detachable Desk Arm,
Swingaway Footrest



Detachable Desk Arm,
Elevating Legrest

Advantage LX		Dimensions and Weights									
		Seat		Back	Overall				Wt. *	Max Wt.	
											
Nylon Upholstery											
3H020100	Nylon, Fixed Full Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H020110	Nylon, Fixed Full Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H020120	Nylon, Detachable Desk Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H020130	Nylon, Detachable Desk Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H020140	Nylon, Detachable Full Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H020200	Nylon, Fixed Full Arm, Swingaway Footrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H020220	Nylon, Detachable Desk Arm, Swingaway Footrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H020230	Nylon, Detachable Desk Arm, Elevating Legrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H020320	Nylon, Detachable Desk Arm, Swingaway Footrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	300 lb
3H020330	Nylon, Detachable Desk Arm, Elevating Legrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	300 lb
3H020340	Nylon, Detachable Full Arm, Swingaway Footrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	300 lb
3H020350	Nylon, Detachable Full Arm, Elevating Legrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	300 lb
Vinyl Upholstery											
3H021100	Vinyl, Fixed Full Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H021110	Vinyl, Fixed Full Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H021120	Vinyl, Detachable Desk Arm, Swingaway Footrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H021130	Vinyl, Detachable Desk Arm, Elevating Legrest	18"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H021220	Vinyl, Detachable Desk Arm, Swingaway Footrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H021230	Vinyl, Detachable Desk Arm, Elevating Legrest	16"	16"	19 1/2"	16 1/2"	32"	36"	26 1/2"	11 1/2"	40 lb	300 lb
3H021320	Vinyl, Detachable Desk Arm, Swingaway Footrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	300 lb
3H021330	Vinyl, Detachable Desk Arm, Elevating Legrest	20"	16"	19 1/2"	16 1/2"	32"	36"	30 1/2"	11 1/2"	40.2 lb	300 lb

*Without front rigging

ADVANTAGE LX – Common Accessories & Parts		HPCPS Code
90763064	Heel loop, pair	E0951
90763062	Wheel lock extension, pair	E0961
B10320F	O2-IV holder	E2208
B10325	O2 cylinder holder	E2208
EJO2BAG07	O2 holder bag	E2208
90763565	Anti-tippers, pair detachable-arm sideframe	E0971
90763566	Anti-tippers, pair fixed-arm sideframe	E0971
92523850	Anti-fold bar (chromed) lower frame	
9B761770	Anti-theft cap for legrest for Traveler chairs with 1.5" spacing	