

Stainless Steel 13 Quart (12.3L) Kick Bucket

Model SS8350

The UMF Medical SS8350 Kick Bucket features all-welded construction, and a non-tip mobile base. A heavy duty rubber bumper protects walls and equipment and also deadens sound. Rubber cushions inside the bucket support ring prevent the pail from rattling.

The UMF Medical SS8350 includes a 13 quart (12.3L) stainless steel pail mounted on a stainless steel mobile frame with four 2" (5.1cm) ball bearing swivel casters.

Comes with a 3-year limited warranty.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › Heavy-duty stainless steel mobile base
- › 13 quart (12.3L) stainless steel pail
- › Rubber cushions inside bucket support ring
- › Rubber bumper on outside of mobile frame
- › 2" (5.1cm) ball bearing swivel casters
- › 3-year limited warranty

Specifications

Height: 14" (35.6cm)
Base Diameter: 14.5" (36.8cm)
Pail capacity: 13 Qt. (12.3L)