

Space-Saving Stainless Steel Foot Stool

Model SS8376

All UMF Medical foot stools boast high-quality, all-welded construction 16 gauge steel. Includes a non-slip corrugated rubber step tread, and no-tip offset stool legs with non-slip rubber tips. All foot stools come with a 3-year limited warranty and a non-corrosive option.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › High-quality, welded stainless steel construction
- › 1.25" x .125" (3.2cm x .3cm) stainless steel band below step platform for rigidity
- › Non-slip rubber feet on stool legs
- › 16 gauge, 1" (2.54cm) outside diameter stainless steel tubing
- › Non-slip rubber step treads
- › 3-year limited warranty

Specifications

Width: 14" (35.6cm)
Height: 7.75" (19.7cm)
Depth: 10" (25.4cm)