

Stainless Steel Stool Transport Cart

Model SS8381

The UMF Medical Stainless Steel Stool Transport Cart boasts high-quality, all-welded construction 16 gauge steel. The convenient push handle makes it easy to transport multiple SS8380 Stackable Stainless Steel Foot Stools at one time.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › High-quality, welded stainless steel construction
- › 16 gauge, 1" (2.54cm) outside diameter stainless steel tubing
- › 1.25" x .125" (3.2cm x .3cm) stainless steel band below step platform for rigidity
- › 2" (5.1cm) swivel casters
- › Convenient push handle
- › 3-year limited warranty

Specifications

Width: 20.5" (52.1cm)
Height (platform): 8.5" (21.6cm)
Height (push handle): 36" (91.4cm)
Depth: 15" (38.1cm)