

new!

gen3

surgical case carts



CDS-333-B

Overall outside dimensions:

29½"W x 28"D x 56"H

Usable Interior Space

20"W x 24¼"D x 44¼"H


One shelf configuration recommended

case carts

surgical

features

- Stainless steel construction.
- Integrated body and floor increases overall cart rigidity.
- High temp, impact resistant bumper mounted directly to cart body reducing the overall footprint.
- Cool-Touch polymer handle.
- 6" diameter stainless steel casters.
- Biohazard Indicator Wheel standard.
- CDS-342 and CDS-345 feature two direction lock casters.
- Roll-out shelves with two position latch mechanism safely locks the shelf in or out.
- Cart is designed to quickly drain during cart washing to minimize residual water.



Cool-Touch Handle



Biohazard Indicator Wheel
CartStat



Modular Bumper



Waterfall Wire Shelf



re•volve™
caster upgrade



CDS-335 shown with optional shelves



CDS-356-MS shown with optional shelves

CDS-335-B

Overall outside dimensions:
35"W x 32³/₄"D x 56"H

Usable Interior Space
25¹/₂"W x 29"D x 44.1¹/₄"H

Two shelf configuration recommended

CDS-342-B

Overall outside dimensions:
46³/₄"W x 28"D x 39"H

Usable Interior Space
37¹/₂"W x 24¹/₄"D x 27"H

One shelf configuration recommended

CDS-345-B

Overall outside dimensions:
46³/₄"W x 28"D x 56"H

Usable Interior Space
37¹/₂"W x 24¹/₄"D x 44¹/₄"H

Two shelf configuration recommended

CDS-356-MS-B

Overall outside dimensions:
29¹/₂"W x 28"D x 39"H

Usable Interior Space
20"W x 24¹/₄"D x 27¹/₄"H

One shelf configuration recommended

Shelving Configurations

CDS-333-SS Stainless steel solid shelf for CDS-333, CDS-356-LD and CDS-356-MS.

CDS-333-SSRO ... Stainless steel roll-out solid shelf for CDS-333, CDS-356-LD and CDS-356-MS.

CDS-333-WS Stainless steel wire shelf for CDS-333, CDS-356-LD and CDS-356-MS.

CDS-333-WSRO .. Stainless steel roll-out wire shelf for CDS-333, CDS-356-LD and CDS-356-MS.

CDS-335-SS Stainless steel solid shelf for CDS-335.

CDS-335-SSRO ... Stainless steel roll-out solid shelf for CDS-335.

CDS-335-WS Stainless steel wire shelf for CDS-335.

CDS-335-WSRO .. Stainless steel roll-out wire shelf for CDS-335.

CDS-342-SS Stainless steel solid shelf for CDS-342 and CDS-345.

CDS-342-SSRO ... Stainless steel roll-out solid shelf for CDS-342 and CDS-345.

CDS-342-WS Stainless steel wire shelf for CDS-342 and CDS-345.

CDS-342-WSRO .. Stainless steel roll-out wire shelf for CDS-342 and CDS-345.

Custom CDS carts engineered to customer specifications available upon request.