



ALEX

- Speech recognition with artificial intelligence (AI) voice responses
- · Scenario-driven, instructor-controlled
- Model-driven physiology with pharmaceutical treatment simulations*
- Virtual AED/Pacing and 12 lead ECG*
- Digital patient and TeamSpace*
- 5-year warranty*
- *Based on service level.

AXEL

- · Scenario-driven, instructor-controlled
- Facilitate patient learner communication via internal microphone and speakers in real-time
- Model-driven physiology with pharmaceutical treatment simulations
- Virtual AED/Pacing and 12 lead ECG
- 5-year warranty

Both ALEX and AXEL are available in light and medium skin tones.



BOTH are ideal untethered/stand-alone simulators to train nurses, medical students and patient care personnel anywhere (no control center required) to:

- Build critical thinking, decision-making & priority-setting skills
- Perform patient physical assessment
- Evaluate vitals and waveforms with simulated patient monitor
- Create automated performance assessment of learning objectives
- Maintain treatment plan with customizable labs and medical orders
- Practice procedures such as I.V. administration, catheterization and airway management
- Assess advanced life support performance with CPR Metrics
- Live stream and/or record patient perspective during simulated events




PC 101005 - 04/21



Nasco Healthcare Inc.

16 Simulaids Drive, Saugerties, NY, 12477 USA Phone: 1-833-NASCOHC (627-2642) Email: info@nascohealthcare.com www.nascohealthcare.com



Nasco HEALTHCARE	ALEXLITE	ALEX PLUS	ALEX PRO	AXEL	
Light Skin Tone	101-7130	101-7140	101-7150	101-7210	_
Medium Skin Tone	101-7130M	101-7140M	101-7150M	101-7210M	
Adult/Adolescent Full Body Patient Simulator	\$	\$\$	\$\$\$	\$\$	
Speech Recognition (NLP)	Up to 250 questions per day	Up to 1,000 questions per day	Unlimited	Manual microphone/speaker allow and student to speak to each	s facilitor other.
Intubation - Oral, Nasal and Digital	Yes	Yes	Yes	Yes	
Palpable Pulses - Bilateral	Yes	Yes	Yes	Yes	
Carotid	Yes	Yes	Yes	Yes	
Brachial	Yes	Yes	Yes	Yes	
Radial	Yes	Yes	Yes	Yes	
Dorsalis Pedis	Yes	Yes	Yes	Yes	
Bilateral Blood Pressure (Manual & 2-step BP) with Smartscope & Smartcuff	Yes	Yes	Yes	Yes	
Anterior Heart and Bowel Sounds with Smartscope	Yes	Yes	Yes	Yes	
Anterior/Posterior Lung Sounds with Smartscope	Yes	Yes	Yes	Yes	
Spontaneous Breathing - Synchronized with Respiratory Rate	Yes	Yes	Yes	Yes	
BVM Ventilation	Yes	Yes	Yes	Yes	
Patient Monitor with Vitals And Waveforms	Yes	Yes	Yes	Yes	
HR, RR, ETCO2, SPO2, TEMP, BP, ECG & MAP	Yes	Yes	Yes	Yes	
/irtual AED/Pacing & 12 Lead ECG	Not Available	Yes	Yes	Yes	
CPR with Metrics Charting Depth, Rate & Ventilations	LOG	LIVE	LIVE	LIVE	
Real-Time Remote Simulation Control	Yes	Yes	Yes	Yes	
Scenarios Included with Customizable Patients	Yes	Yes	Yes	Yes	simble (Elles (20 km)
Measurable Objectives with Automated Assessment	Yes	Yes	Yes		eight: 45 lbs. (20 kg) eight: 5 ft. 1 in. (1.55 m)
Patient Orders	Yes	Yes	Yes		Shipping Weight: 67 lbs. (Shipping Dimensions:
Model-Driven Physiology	Not Available	Not Available	Yes	Yes Sh	
Drug Treatments	Not Available	Not Available	Yes	Yes 41	" H x 19.5" W x 15" D
Built-In Microphone and Speakers	Yes	Yes	Yes	Yes	
Events, Objectives & Orders Log Entry Record	Yes	Yes	Yes	Yes	
Communication Transcription Log Entry Record	Yes	Yes	Yes	Not Available	
ris Cam - Live Stream	Yes	Yes	Yes	Yes	
ris Cam - Recording	Not Available	1,000 hours	Unlimited	1,000 hours	
Simulation Log Archive	60 days	1 year	5 years	1 year	
eamspace	Not Available	Not Available	Yes	Not Available	
Cloud Connected Control	Yes	Yes	Yes	Yes	
/irtual Patient	Not Available	Not Available	Yes	Not Available	
Warranty	1 year	5 years	5 years	5 years	
V Arm - Optional	Optional	Optional	Optional	Optional	
Female Conversion Kit - Optional	Optional	Optional	Optional	Optional	

ACCESSORIES

ALEX/AXEL Female Accessory Kit Light Skin Tone ALEX/AXEL Female Accessory Kit Medium Skin Tone ALEX/AXEL Optional IV Arm Light Skin Tone - Right 101-7160 101-7160M 101-8029

101-8029M ALEX/AXEL Optional IV Arm Medium Skin Tone - Right 101-8030UPPLUS ALEX Upgrade Lite To Plus 101-8030UPPRO ALEX Upgrade Lite To Pro

101-8040UPPRO ALEX Upgrade Plus To Pro 101-8040EXS ALEX Plus Extended Yearly Service Fee 101-8050EXS ALEX Pro Extended Yearly Service Fee